



First Medical Source

Reducing Cost of Healthcare Worldwide



FMS Dual Post-Operative Local Anesthesia Kit
High Quality
Easy to use
Not made with natural rubber Latex or DEHP
FDA compliant, and Market tested
Responsibly Priced

Indication for use

The FMS Local Anesthesia system is ideal for regional or local Anesthesia for control of acute post operative pain. The FMS pump delivers a continuous fixed rate infusion of non-narcotic local anesthetic to the surgical site through a catheter which is inserted at the site during surgery. The FMS pump is safe, easy to fill, simple to use and lightweight, so that patients can keep their mobility during their recovery.

First Medical Source is proud to offer a safe and less expensive alternative to existing Local Anesthesia kits in the market. Our expertise and knowledge of the supply chain enable us to control cost, quality, and delivery to support your need. Please review the available models of the FMS Local Anesthesia pumps and catheter kits described below and then contact First Medical Source to order or to obtain more information.

Regional Anesthesia Kit Catalog



To place an order, send an email to Sales@firstmedicalsourc.com or call (949) 637-5222. For additional information, email us at Info@firstmedicalsourc.com

To report a problem or file a complaint, email us at Support@firstmedicalsourc.com



First Medical Source

Reducing Cost of Healthcare Worldwide

www.firstmedicalsourc.com



Manufactured for and distributed by: First Medical Source, LLC,
Laguna Niguel, CA 92677

© 2008-2014 by First Medical Source, LLC, All rights reserved.
PMK-1102H H 01/2014

Order your selected FMS™ Pump and Catheter Kits

Catheter Kit Options



Catheter Kit Description	Kit Order Code
2.5" (15 holes) Fixed-Hub FMS Conduction Catheter™	
Single Site with a 4.5" Tear-Away Introducer	SC25
Dual Site with a Y-Connector & two ea 4.5" Tear-Away Introducers	YC25
5" (30 holes) Fixed-Hub FMS Conduction Catheter™	
Single Site with a 4.5" Tear-Away Introducer	SC50
Dual Site with a Y-Connector and two ea 4.5" Tear-Away Introducers	YC50
10" (60 holes) Fixed-Hub FMS Conduction Catheter™	
Single Site with a Tear-Away Introducer - sold separately	SC100

Each catheter pack includes:

- FMS Catheter with female luer cap
- Introducer
- 60 cc filling syringe
- 5 cc priming/flushing syringe
- Transparent dressing
- Adhesive strip cards
- Bifurcated set (dual packs only)

Rx CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Introducers and Pump Options

Product Description	Order Code
FMS Single Pack Sterile Tear-Away Introducers	
	
8" 17 G Introducer	INT8
12" 17 G Introducer	INT12
Infusion Pumps	
	
250 ml volume at 5.2 ml per hr continuous (2 days)	C0052L
250 ml volume at 2.1 ml per her continuous (5 days)	C0021L
100 ml volume at 2 ml per hr continuous (2 days)	100D2
Other models can be supplied by custom order.	TBD